



**MODEL**  
**F150A**  
Dental Intra-Oral Sensor  
1.5 Size

**MODEL**  
**F350A**  
Dental Intra-Oral Sensor  
2.0 Size

**PRODUCT FEATURES**

- Compact head size
- Ergonomic design that minimizes foreign body feeling in the mouth
- High-definition image with 14bits gray scale and best actual resolution
- 3 types of preset image mode for user convenience
- Printing mode with reporting function
- Easy to install software

**PRODUCT SPECIFICATIONS**

Model	F350A	F150A
Dimensions (W x L)	33.7 × 43.4 mm	27 × 37 mm
Scintillator Technology(CsI)	FOS	FOS
Active Area(CMOS)	26 × 34 mm	21.6 × 30 mm
Theoretical Resolution	25 lp/mm	
Pixel Size	20um	
Gray Levels	16,384(@14bits)	
Cable Length	2.5m	
Thickness	5.4mm	
Connection	Standard USB 2.0	

**CERTIFICATION STATUS**

**ISO 13485**

**GMP**

**FDA510 (K)**

**MFDS**

**INTRA\_ORAL SENSOR SOFTWARE (CONSOLE VIEWER)**

