



**MODEL**  
**P70A**  
Dental Portable X-ray


**MODEL**  
**PDXT**  
70kV/2mA  
CNT Digital X-ray Tube

**PRODUCT FEATURES**

- CNT digital x-ray tube
- Perfect repeatability
- No stray X-rays (Totally low dose)
- No oil leakage
- Ergonomic and intuitive interfacing

**CERTIFICATION STATUS**

• **ISO 13485** • **GMP** • **FDA510(K)** • **MFDS**



**PRODUCT SPECIFICATIONS**

Tube Voltage	70kV	Exposure Time	0.05~1.0s	Total Weight	1.6kg
Tube Current	2mA	Lifespan	200,000 shot	S.S.D.	200mm
Focal Spot (nominal)	0.4	Maximum Duty Cycle	1:30		